

# Generation and Analysis of Biolog Data for Metabolic Modeling



Heather Danhof, Firas Midani, James Collins

Laboratory of Robert Britton

Baylor College of Medicine, Houston

# Overview:

- Introduction to what the Biolog system is
- Rationale for protocol modifications
- Example data and talking points
- Data workflow and analysis.






# Phenotype MicroArray Technology

Patented 96-well  
array system

Standardized  
comparisons  
across strains

Medium through  
put



  
**BiOLOG**

# How the Biolog system works

- Utilizes proprietary reporter system of NADH production to reduce a tetrazolium dye and produce a color change (purple).
- The regents are designed to slow respiration in “weakly” positive cells and reduce any false positive signals.
- Uses a specialized OmniLog instrument to read and analyze data



- NO Gram stain
- NO pre-tests
- NO follow-on tests
- ONE panel for both GN & GP bacteria
- ONE minute set-up
- OVER 1350 species coverage

<https://www.toshindia.com/products/bacteria-yeast-and-fungi-identification-system>

# Using Biolog Phenotype Arrays

## **Benefits**

- High throughput as marketed
- Convenient pre-arrayed compounds
- Standardized

## **Limitations**

- Proprietary
  - No details on concentrations
  - No ability to change media components
- Reagents and dyes are very oxygen sensitive
- Special equipment is costly

# Modifications to the Standard Biolog Protocol

- OD<sub>620</sub> Absorbance to monitor growth directly.
  - Uses a standard plate reader inside the chamber.
    - Programmable readings
    - Control over chamber environment
  - Solves the issue of the oxygen sensitive dyes in the anaerobic chamber.
- Use completely defined minimal growth medium that has been optimized for *C. difficile*.



# Experimental setup

Strains are grown overnight in BHIS

Diluted 1:10 in CDMM

Aliquoted in Biolog plates (100 uL/well)

Incubated in anaerobic plate reader, 37 degrees


Shakes 30 seconds before reading every 10  
minutes for ~16 hours

Each strain is done in duplicate, sometimes more\*





# Biolog data


Each plate has unique growth curves




# Plate to plate variability




## Max OD600



## Fold Change



# Determining Growth and fold change



$$\text{Fold change} = \frac{\text{Max Growth Xpt. Carbon}}{\text{Max Growth Neg. control}}$$



Press ESC to cancel Copy

# Reproducibility

Xpt Compound

L-Arabanose

N-Acetyl-D-Glucosamine

D-Saccharic Acid

Succinic Acid

D-Galactose

L-Aspartic Acid

L-Proline

D-Alanine

D-Trehalose

D-Manose

Dulcitol

D-Serine

D-Sorbitol

Glycerol

L-Fucose

D-Glucuronic Acid

D-Gluconic Acid

D,L-alpha-GlycerolPhosphate

D-Xylose

L-Lactic Acid

Formic Acid

D-Mannitol

L-Glutamic Acid

D-Glucose-6-Phosphate

D-Galactonic Acid-gamma-Lactone

D,L-Malic Acid

D-Ribose

Tween 20

L-Rhamnose


D-Fructose

Acetic Acid

alpha-D-Glucose

	Repeated	Repeated
Sym	Sym	Sym
PRB807	PRB807	PRB807
RT002	RT002	RT002
S0065	S0065	S0065
1.2	0.9	0.9
3.1	3.2	2.8
1.1	1.1	0.9
1.4	1.4	1.2
1.0	1.1	0.9
1.0	1.1	0.9
1.6	1.0	0.8
1.0	1.0	0.8
1.3	2.3	1.9
3.2	2.4	2.0
1.0	1.1	0.9
1.2	1.2	1.0
3.6	2.3	1.9
1.0	1.1	1.0
1.0	1.1	0.9
0.9	1.1	1.0
1.1	1.1	0.9
1.0	1.1	0.9
1.2	1.2	1.3
1.0	1.1	0.9
1.0	1.1	0.9
3.3	2.8	2.2
1.1	1.1	0.9
1.1	1.1	0.9
1.1	1.1	1.0
1.1	1.1	0.9
1.1	1.2	1.3
1.0	1.0	0.8
1.0	1.0	0.9
3.4	3.3	2.7

D-Sorbitol



***Comparative analysis:***

***Comparative analysis:***

- (i) Growth on same carbon substrate across different genotypes.*



***Comparative analysis:***

- (i) Growth on same carbon substrate across different genotypes.*
- (ii) Growth of isolates from same genotype on different carbon substrates.*

***Comparative analysis:***

- (i) Growth on same carbon substrate across different genotypes.*
- (ii) Growth of isolates from same genotype on different carbon substrates.*

***General questions:***

- 1. Does each isolate have a unique carbon substrate utilization profile?*

## ***Comparative analysis:***

- (i) Growth on same carbon substrate across different genotypes.*
- (ii) Growth of isolates from same genotype on different carbon substrates.*

## ***General questions:***

- 1. Does each isolate have a unique carbon substrate utilization profile?*
- 2. Do isolates of the same genotype grow on the same carbon substrates?*

## ***Comparative analysis:***

- (i) Growth on same carbon substrate across different genotypes.*
- (ii) Growth of isolates from same genotype on different carbon substrates.*

## ***General questions:***

- 1. Does each isolate have a unique carbon substrate utilization profile?*
- 2. Do isolates of the same genotype grow on the same carbon substrates?*
- 3. Is there strong correspondence between genotype and phenotype?*




## **Comparative analysis:**

- (i) Growth on same carbon substrate across different genotypes.*
- (ii) Growth of isolates from same genotype on different carbon substrates.*


## **General questions:**

- 1. Does each isolate have a unique carbon substrate utilization profile?*
- 2. Do isolates of the same genotype grow on the same carbon substrates?*
- 3. Is there strong correspondence between genotype and phenotype?*
- 4. Which substrates may be driving adaptation in C. difficile populations?*

# Overview of data processing




## Overview of data processing



*Simple implementation is possible with a scientific computing language (e.g. Python or R) and/or a spreadsheet (e.g. Microsoft Excel).*

# Overview of data processing






# *Biolog Phenotype Microarrays capture utilization of carbon substrates.*

PRB954 PM1


Fold change > 1.5 or Fold change < 0.5




# Manual and automated quality checks are necessary!

CD801 PM1

Fold change > 1.5 or Fold change < 0.5




# Overview of data processing




# *Biolog Phenotype Microarrays capture utilization of carbon substrates.*

PRB954 PM1


Fold change > 1.5 or Fold change < 0.5



# *Clinical isolates exhibit varying growth dynamics on different carbon substrates*





### *Clinical isolates exhibit varying growth dynamics on different carbon substrates*






*Final time points do not reflect the true relative fitness of a carbon substrate.*




*Final time points do not reflect the true relative fitness of a carbon substrate.*




*There is clear trade-off between assay  
throughput and biological insight!*

### *Clinical isolates exhibit varying growth dynamics on different carbon substrates*




# *Varying growth dynamics occur also on the same carbon substrate*


## Growth on Trehalose




# Overview of data processing





*Growth curve dynamics are more informative than final OD measurement.*



*A non-parametric approach for growth curve fitting may be useful.*




*A non-parametric approach for growth curve fitting may be useful.*




Tonner et al. *Genome Biology*. 2017.



# Overview of data processing




*Baseline growth of different genotypes can vary and bias downstream analysis.*



Subset of carbon substrates that are utilized by at least two isolates.

*Baseline growth of different genotypes can vary and bias downstream analysis.*




*Instead of a fold change that can range between 0 and  $\infty$ , we can transform and scale the data for example into z-scores.*

$$Z_c = \frac{(F_c - \mu_g)}{\sigma_g}$$

*where  $F_c$  is the fold change on a specific carbon, while  $\mu_g$  and  $\sigma_g$  are the mean and standard deviation, respectively, of fold changes for a specific genotype.*


Subset of carbon substrates that are utilized by at least two isolates.

# How well does phenotype correspond to genotype?




Subset of carbon substrates that are utilized by at least two isolates.

# *How well does phenotype correspond to genotype?*





*Which substrates may explain most of the variation between the different isolates?*



Feedback is welcome on how improvements to data workflow can better address questions & objectives.

Many ways to analyze the data,  
but outcome strongly depends on  
the input and data manipulation.



*Feedback is welcome on how improvements to data workflow can better address questions & objectives.*

**Comparative analysis:**

- (i) Growth on same carbon substrate across different genotypes.*
- (ii) Growth of isolates from same genotype on different carbon substrates.*

**General questions:**

- 1. Does each isolate have a unique carbon substrate utilization profile?*
- 2. Do isolates of the same genotype grow on the same carbon substrates?*
- 3. Is there strong correspondence between genotype and phenotype?*
- 4. Which substrates may be driving adaptation in *C. difficile* populations?*

